
Gebrauchsanleitung
Instructions for use
Instructions pour l'usage
Istruzioni per l'uso
Modo de empleo

Ε 90

infrarot

E 90 Infrarot-Mono-Empfänger

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E 90 infrared mono receiver

The wireless infrared mono receiver E 90 receives music and speech radiated by an appropriate infrared transmitter.

Types of application:

- TV and hi-fi sound reproduction
- as a delegate headphone, for interpreter systems
- in theatres or during conferences as a hearing aid

Features:

- hi-fi quality
- extremely light in weight
- easy to handle
- economic accu operation causing no harm to the environment
- recharging of BA 90 accu in the charging port of the transmitter or by means of the L 90 charging unit

Technical data:

Receiving frequency	95 kHz \pm 25 kHz, FM
Frequency range	30-18,000 Hz
S/N ratio	65 dB(A)
THD at 1 kHz	< 1 %
Max. sound pressure level at 1 kHz	approx. 115 dB
Power supply	BA 90 accu (Item No. 03261)
Operating time	approx. 11 hrs
Accu recharging time	approx. 18 hrs
Weight	approx. 40 g (including accu plug)

Fault finding:

Hissing noise only	Transmitter is not set into operation or it is defective
Volume turned to max. but no signal available	Signal is missing, check connection to transmitter! Volume can be adjusted separately at headphone socket? See manual of TV/hi-fi components!
Sound is disturbed	Accu is completely discharged, please recharge it immediately!
Sound is slightly noisy	The transmitter is not positioned in an optimum way, reception is obstructed.

Before using it for the first time, the accu has to be charged for 24 hrs!
Please switch off the device after use!

Ear cushions (removable) have to be cleaned with soap-suds from time to time!

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